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### EN Instructions for training the pelvic floor with Joyballs secret

Congratulations, you have decided to purchase the innovative Joyballs secret pelvic floor training balls, which are distinguished by not only a first-class design, but also the following features:

- patented retrieval tab which is worn discreetly inside the body (Joyballs secret)
- whisper-quiet thanks to patented technology (Joyballs secret)
- trampoline effect\* for gentle stimulation and more effective training of the pelvic floor muscles

\*The trampoline effect is an internally concealed innovation in which the structure of the balls' interior walls was constructed in such a way as to create random changes in the direction the balls move, which is considerably better for pelvic floor training than a continuously one-directional movement. That is because such random movement provides

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permanently new stimulation impulses (similar to electric stimulation therapy). We call this intelligent effect the "trampoline effect" – for more effective training, which in practice translates as "more effect for less movement".

Your Joyballs secret will provide you with reliable and discreet support in the effective prophylactic training of and therapy for your pelvic floor muscles (musculus pubococcygeus). Joyballs secret are a branded product which strengthens your pelvic floor muscles simply, effectively, and discreetly. Training with Joyballs is recommended by gynaecologists and midwives.

The pelvic floor muscles serve an important function within the female body. They hold the pelvic organs in their proper, functional position. Pregnancy, vaginal births, and other bodily stresses can overstretch the pelvic floor muscles, and result in the malposition of the pelvic organs, e.g. dysfunction of the bladder, prolapse of the uterus (descensus uteri), etc.

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**EN** A prolapse of the vagina (genital prolapse) is an extreme example of the disorders which can be caused by pregnancy and giving birth, injuries sustained while giving birth, great physical stresses, or congenital weakness of the connective tissue. The cause of the prolapse of either the vagina or the womb is a weakness of the connective tissue in the pelvic and/or an over-stretching and weakening of the pelvic floor muscles. This prolapsing varies in its level of severity. What occurs is a complete or partial protrusion of parts of the vagina and/or womb.

Despite changes in social attitudes and increased awareness concerning health and prevention, certain areas of the body such as the genital area are still associated with strong taboos and shame. It is therefore not surprising that the pelvic floor still languishes in relative obscurity.

One reason for this may lie in socially conditioned rules and patterns of behaviour. Even in the 21st century, people are still lacking in knowledge about their own intimate areas. Another reason seems to be that the pelvic floor muscles are not visible, i.e.

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are "hidden" within the body. It's not that easy to feel the right muscles and train them in a targeted manner.

The pelvic floor often receives attention only once the weakened muscles of the pelvic floor have caused noticeable and oftentimes worrying consequences. One of the dreaded consequences is urinary incontinence. A healthy, strong pelvic floor is a prerequisite to many processes within our bodies. It carries and supports the abdominal organs and allows for the voluntary activation of the sphincter muscles. It makes a decisive contribution to a proper sense of wellbeing and to good posture.

Targeted training through the use of Joyballs secret can restore strength the pelvic floor muscles, correcting or preventing the dysfunction of the bladder and other pelvic organs.

A strong pelvic floor reduces the risk of prolapses, and can increase both partners' pleasure and the ability to climax during sexual intercourse.

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**EN** With our guidebook, we aim to provide you with an insight into this muscle group, increase your awareness, and explain the links between your pelvic floor and your general wellbeing. We will also show you practical exercises and tips for strengthening and looking after your pelvic floor.

This guidebook is not intended as a handbook for self-treatment. Rather, it's intended to boost your awareness and self-initiative, and to encourage you to become pro-active. It is vital to obtain a thorough clinical assessment of the causes .e.g. leading to incontinence, as well as a clear diagnosis.

We recommend that you consult a specialist and obtain from him or her extensive instructions on pelvic floor training. Physiotherapists and midwives who have undergone additional training are particularly appropriate points of contact, including regarding post-natal therapy. Pelvic floor training always requires individual motivation, regularity, and perseverance. Take a little time to train, and your body will thank you for it.

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Exercising with the Joyballs secret is simple, and can yield successful results after just a few weeks of consistent training. A prerequisite for this is that the exercises are carried out regularly and carefully.

Joyballs help you increase sensation in your pelvic floor by enabling you to actively train it. Even when gymnastic training of the pelvic floor alone has failed owing to an inability to feel the pelvic floor muscles, the combined training with Joyballs can lead to the desired therapeutic success.

Exercising with the Joyballs is especially suitable for the following indications:

- strengthening the musculature following child-birth
- correcting urinary incontinence and/or double incontinence in women caused by a weak pelvic floor
- preventing a prolapse of the uterus (descensus uteri) and other pelvic organs
- after surgery to the lesser pelvis

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### EN Structure of the pelvic floor

The various muscles of the pelvic floor are often not only unfamiliar and unfeelt, but also weak. In order to train the pelvic floor, it is useful to know a few basic facts about its structure and function.

What does the pelvic floor look like?

The pelvic floor consists of three layers of muscle. The musculature of the pelvic floor is a complex, multi-layered structure, with layers running horizontally, vertically, and diagonally. It has an interactive relationship with the abdominal organs (bladder, intestines (rectum), womb), as well as affecting posture and stability, breathing, and, not least, sexuality.

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Outer layer – the eight shape:

The outer layer surrounds and strengthens the pelvic floor's orifices around the anus, vagina, and urethra, and looks like an elongated figure eight. The muscle fibres of this layer originate at the front (pubic bone) and insert at the rear (tail bone, or coccyx). The condition or strength of this layer is of primary importance to continence, as it supports the outer sphincter muscles.

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
### EN Inner layer – the fan:

The second, inner layer covers a greater area, from the lower back and tail bone to the sides of the pelvis. The muscles form a V-shape. At the centre of the V, the opening of the anus, urethra, and vagina remain free. The levator ani muscle is located down in the bony pelvic brim (like a bowl shaped out of muscle).

### Middle layer – the triangle:

The middle layer forms the front part of the pelvic floor and is positioned diagonally to the two other layers, so as to stabilise the opening in the middle. It links the ridges of the ischial tuberosity at the level of the perineum, i.e. the area between the anus and the sexual organs.

The muscle groups of the pelvic floor are cross-linked with one another through muscle fibres and networks of nerves. This leads to our contracting and relaxing all the muscles simultaneously the majority of the time.

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The pelvic floor muscles play a defining role in the posture of the spine. Sufficient tensile strength in the pelvic floor prevents or at least compensates for poor posture in the spine, pelvis, sacrum, and head.

It is important to know that the pelvic floor represents a system that functions in unison. The pelvic floor is like the elastic membrane of a trampoline. It consists of a layer of muscle that runs through the rectum, vagina, and urethra.

There are also retaining bands along the edge, which oscillate like springs.

The muscle layer (membrane) and the elastic retaining bands are anchored to the pelvis at the front, middle, and rear. The womb and bladder are positioned on top of the pelvic floor.

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**EN** ... and how it works:

The vagina acts like an elastic membrane.

Contracting the front area of the vagina seals the urethra; contracting the rear area seals the bladder. This allows you to retain urine.

If the connective tissue of the vagina or the retaining bands weakens, it becomes impossible to seal the urethra. So-called stress incontinence occurs, leading to involuntary urination (e.g. when coughing).

This same weakness of the connective tissue is often the cause of instability in the bladder. This causes the lower portion of the bladder to be prematurely stimulated, even when not properly full.

The result is a constant urge to urinate, day and night (urge incontinence).

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## The pelvic floor performs three major functions:

### Contraction

This plays a role in continence, as described previously. Contraction of the pelvic floor allows a person to voluntarily pass urine and stools.

### Relaxation

When emptying the bladder or rectum, the pelvic floor relaxes so that the sphincters of the urethra and anus can open.

### Reflexive resistance

When coughing, laughing, or carrying heavy loads, the abdominal area, and hence the pelvic floor is contracted. The pelvic floor contracts even further during these activities to prevent the passing of any urine or stool.

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## EN The ageing pelvic floor

The pelvic floor ages the same way as the rest of the human body, the organs shrink and the retaining bands lose their elasticity. The pelvic floor becomes less elastic and its retaining strength weakens.

This makes pelvic floor training whilst still younger all the more important. A slight weakening of the pelvic floor muscles is as yet not great cause for concern. However, if the ageing of the pelvic floor is already advanced, the previously described complaints can arise. Urinary and faecal incontinence can occur, and the womb may prolapse, which can have extremely unpleasant consequences.

### Risk factors

There are a few risk factors that make a weakening of the pelvic floor's supportive properties more likely:

- an untrained pelvic floor: an untrained pelvic floor provides little resistance to the ageing process or the previously described risk factors; it easily loses its stability.
- multiple births: the birth of each child places extreme stress on the pelvic floor; if it is not sub-

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sequently trained, its supportive properties greatly diminish in the course of time.

- frequent lifting and carrying of heavy objects: this also places a disproportionate amount of stress on the pelvic floor.
- hormones: when a woman's oestrogen production levels are reduced with the onset of the menopause, this also has a negative effect on the pelvic floor muscles.
- being overweight: the heavier a person is, the harder the pelvic floor has to work to be able to fulfil its functions.
- age: the older a person is, the more the muscle weakens.

### The taboo topic of bladder weakness

The problems experienced do not necessarily have to involve (just) urinary incontinence. One can also develop problems such as pressure in the vagina, or deep-seated pain in the sacrum. Many people then suspect they have arthritis or a slipped disc, yet no doctor can find the problem and prescribe the necessary treatment.

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**EN** The sexual functions can also be disrupted, from a loss of libido through to real pain which makes everything impossible.

Women often suffer from recurring bladder infections, and are occasionally also unable to properly control their defecation. This results in either faecal incontinence of varying degrees of severity or constipation with the formation of haemorrhoids.

The causes can vary; a difficult birth which causes the tearing of the pelvic floor's retaining bands; congenitally weak connective tissue, and disproportionate stressing of the pelvic floor through excess body weight, through coughing during chronic lung diseases, and through hard physical labour.

What are the Joyballs made of?

For even greater safety and comfort, our easy-to-clean Joyballs are manufactured using SILIKOMED®, a hypoallergenic as well as clinically and dermatologically tested material.

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How do I clean the Joyballs?

It is vital that the Joyballs be thoroughly cleaned before and after each use. Warm water and soap are perfectly sufficient for cleaning your Joyballs.

We recommend the use of clean'n'safe, a disinfectant that has been proven to work excellently with Joyballs.

When should Joyballs not be used?

- when your doctor advises against their use
- during menstruation
- in the event of vaginal inflammation
- in the event of a vaginal infection
- when recovering from any injury to the vagina (e.g. following an episiotomy)
- you should use Joyballs during and following pregnancy only after consulting your gynaecologist or midwife

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### EN How do I train?

#### Pelvic floor training

The American gynaecologist Dr. Arnold H. Kegell (1894–1981) “invented” pelvic floor training. Dr Kegell was a gynaecologist who noticed that many women suffered from temporary urinary incontinence problems following childbirth. He therefore developed a series of exercises to strengthen the pelvic floor muscles. This why exercising the pelvic floor is also called “Kegeling”, in honour of its “inventor”.

Pelvic floor gymnastics are simple, and what’s really handy about them is that many pelvic floor exercises can be carried out anywhere, without anybody realising what you’re doing. However, it’s important to carry out the exercises carefully and correctly.

#### There are two forms of training:

The first type is pure bio-feedback training with vaginal training balls.

Firstly, this allows you to develop an awareness of how to feel your pelvic floor muscles. Secondly, when you first begin with the training and your pel-

vic floor muscles are as yet untrained, the Joyballs work on a “bio-feedback” principle – a type of reflexive reaction from the muscles. Each ball contains a weight, which reacts to the movements of your body. Owing to their own weight (and the ergonomically adapted design of the Joyballs secret), the Joyballs try to slip out – which your pelvic muscles try to prevent. The muscle contractions during a training session function like those performed during strength training with weights. Moreover, the balls inside the Joyballs oscillate against the outer wall with each movement. Thanks to an innovative integrated “trampoline effect” (for which the Joyballs secret are especially optimised), gentle “kinetic vibrations” are set in motion, which both encourage your pelvic floor muscles to contract and relax, and provide you with stimulation.

A single training session should last up to 10 minutes. Training can be performed while either standing still or walking around. Two daily training sessions of up to 10 minutes should be sufficient initially.

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**EN** Once you can hold the pelvic floor training balls in for over 10 minutes without difficulty, you should increase the duration of training and combine it with gymnastic pelvic floor training.

The second type is combined bio-feedback gymnastics.

Here, the benefits of the previously described bio-feedback training is combined with special training exercises for the pelvic floor. Combined bio-feedback gymnastics represent the most effective training method for achieving optimal results.

The position of the pelvic floor muscles

It is important to find the right muscles. Studies have shown that 30% of all women trying to train their pelvic floor muscles are doing so incorrectly. Common errors include contracting the buttocks instead of the pelvic floor muscles, holding one's breath, and pushing downwards rather than squeezing and lifting. Differentiating between the muscles of the buttocks and the pelvic floor: often, and especially when starting out, the buttock mus-

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also are also contracted when training. This doesn't matter, but you should be aware that this is not the same muscle group.

To recognise the difference between the muscles of the buttocks and those of the pelvic floor, sit on a hard chair and lean forward slightly. Now strongly contract all the muscles around your genital area: when you contract the muscles, you will feel a light upward movement of the body, followed by a downward movement when you relax: this shows that you are using the buttock muscles. Now try contracting the muscles in your genital area without causing the upward and downward motion just described. With time, you will become ever better at differentiating between the different muscle groups.

Many people have little awareness of this area of muscle. Important is initially to find and feel the muscles in your pelvic floor. The great advantage of pelvic floor balls is the improved perception of the correct muscles.

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
**EN** And before you begin, think about the proper way to breathe: exhale when contracting the pelvic floor! Inhale when relaxing the pelvic floor!

The following test should be helpful for this. When urinating, try to stop the flow mid-stream. Feel which muscles you are employing for this. This is the muscle group you need to train. Be aware that this is just a test for finding the right muscles. You should not interrupt the flow of urine doing urination.

When you contract the pelvic floor muscles correctly, you will feel pressure and a slight upward motion extending from underneath the pelvis to the inside of the body, around the vagina and anus. Other body parts should not move at this point (e.g. the buttocks or the inner thighs). Continue breathing steadily.

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Inserting pelvic floor training balls in example of the Joyballs  
Now gently insert the Joyballs into the vagina (like a tampon). We recommend the use of a lubricant (e.g. AQUAglide or BIOglide) for easy and pain-free insertion. The Joyballs should sit deeply and almost vertically, with the retrieval tail facing downwards (see image).

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## EN Strengthening the musculature

Once you have learned how to use these muscles properly, train them regularly – ideally, from one time up to a maximum of three times a day. Daily training is recommended in order to ensure lasting training progress. It takes only a few minutes and can be carried out in various positions (e.g. when sitting, standing, or lying down). No one can see what you are doing. Initially, you may only be capable of maintaining the contraction for one to two seconds. Increase this to a maximum of ten seconds. You should rest in between each contraction, with the rest lasting at least as long as the contraction (e.g. five seconds of contracting, five seconds' rest). Initially, you may feel tired after just a few contractions. Gradually increase your performance to ten repetitions. You can also train using a few quick, powerful contractions. Contract the muscles as powerfully as possible before relaxing them again. Repeat this exercise up to ten times.

Make use of these muscles when you need them. In conjunction with these exercises, you should also

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contract your pelvic floor muscles when you need them, e.g. when laughing, coughing, bending over (situations which cause urinary incontinence). For optimal results, you should continue with the training described here for a period of around six months. After two to four months of regular exercising, you should notice an improvement. Once you have achieved the desired effect, you can decrease the frequency of training. However, make sure you use these muscles as described every day when you require them. Research has shown that women who have trained intensively retain the improvements for months and even years. It is never too late to begin pelvic floor training. Women in their seventies and eighties can improve and even eliminate their medical complaints. If you integrate these exercises into your daily routines, they will become as automatic as brushing your teeth. If, despite conscientious training, you still do not notice any improvement, please consult your doctor. He or she can advise you on additional measures and, where appropriate, refer you to a specialised physiotherapist or a gynaecologist/urologist.

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**EN** When you first begin training, it can leave your muscles feeling sore in the area of your pelvis and bowels. This is a sign that you have carried out the exercise correctly, and that previously untrained muscle has been strengthened. – You should therefore be careful not to train for too long initially, instead gradually and continuously increasing the duration of training.

In general you can carry out pelvic floor exercises when seated, standing, or lying down. As described above, you should contract the pelvic floor muscles as hard as possible without using any other muscle group while doing so.

Tension and relaxation – alternation brings success!

It is only the alternation between contraction and relaxation that creates the actual training effect. During each exercise, ensure that you contract the muscles as strongly as possible, and maintain the tension for approximately 6–8 seconds. Breathe normally whilst doing so, and do not tense the sto

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mach muscles. N.B.: After each contraction, rest for at least as twice as long as the contraction phase!

Repeat this exercise up to ten times. Try to contract the muscles for six to eight seconds each time. You should carry out the exercises at least three times each day.

It will be even more effective if you first contract the muscles for eight seconds and then try to contract the musculature even further with three to four rapid contractions.

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## EN Preliminary exercise: Feeling for the pelvic floor

It isn't always easy to find the precise location of the pelvic floor.

You should therefore perform the following exercise:

Sit on the floor and cross your legs. Sit up straight, look straight ahead, and ensure that you are sitting on both your sitting bones (the ischial tuberosity). Place both your hands on your lower abdomen and close your eyes. Then try to push together the sit bones without moving. At the same time, contract your muscles in the way you would if you desperately needed to urinate, and pull your navel inwards and upwards. Maintain the tension for 15 to 20 seconds, then slowly release it again. While wearing pelvic floor training balls contract the muscles around the balls.

Repeat this exercise several times until you can feel the contraction and relaxation of your pelvic floor.

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
## EN Posture exercises

It is important to maintain an upright, pelvis-friendly position in daily life. This exercises will help you to improve your awareness:

Take a hip-width stance, bending your knees only slightly. Straighten your spine and tilt your head upwards from a point originating from the crown. Place one hand on your lower abdomen, the other on your chest. Straighten your sternum and tilt it diagonally upwards. Pull your navel inwards and upwards whilst pulling your shoulders back and down. Stay in this position for a few breaths, then relax.

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## EN Back strengthener

Stand in an upright position, elongate your spine, and maintain a neutral head position. Take a shoulder-width stance, with the knees and toes pointing forward. Raise your arms and position your elbows at shoulder height, with an approximately 90 degree angle (right angle) between the upper arm and the forearm. Pull your shoulders back and down. Activate your pelvic floor, and tilt your torso forwards from the hips, until the angle between your thighs and your torso is around 90 degrees. Slowly stretch your arms out in front of you. Then bend your arms again. Repeat this 25 times before slowly relaxing.

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
## EN Bending exercise

Stand with your legs slightly apart and bent. Prop your hands on your thighs whilst keeping the back straight. Pull your pelvic floor muscles upwards and inwards.

Repeat this exercise at least eight times.

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
## EN Walking exercise

Stand upright and resume tension in your pelvic floor muscles when exhaling. Release it once more when breathing in. Now try to maintain the tension when breathing in as well, and at the same time walk on the spot for a few steps. Maintain the tension for as long as you can whilst avoiding any unnecessary use of the buttock or abdominal muscles. Continue to train your pelvic floor by consciously squeezing together your pelvic floor when walking. Don't stop your breathing rhythm!

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
## EN For "advanced trainees" – the virtual lift

Stand upright with a hip-width stance, your knees unlocked. Now contract your pelvic floor muscles in stages: Tension is increased as in a lift that is going up and stopping at each floor on the way. Try to briefly maintain the tension level "at each floor" and continue to breathe evenly before gradually releasing the tension, i.e. letting the lift go "down to the ground floor". Once the ground floor has been reached, the lift continues to the "basement level" for deep relaxation of the pelvic floor. Once you have become familiar with this exercise, you can increase the height of your virtual "exercise building" from 3 to 4, 5, 6, and more floors. Think about "sucking in" and don't contract your buttock or abdominal muscles!

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
## EN Diagonal crunch

Lie on your back, pull your knees up, and place your feet flat on the floor to begin with. Put a folded-up towel, a blanket, or a small cushion between your knees and slowly lower your knees to the left. Place your hands on the back of your head and let your elbows fall to either side. Squeeze your knees together and activate the tension in your abdominal muscles. Slightly raise your head and shoulders off the floor. From this position, gradually raise your upper body to the fullest extent possible without using any rocking motion. Slowly relax again.

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
## EN Lower abdominal lift

Lie flat on your back. Put a folded-up towel, blanket, or small cushion between your knees. Put your arms to your sides and relax them. Elongate your spine. Gently contract your stomach and pelvic floor muscles and lift your pelvis off the floor. Firmly contract your pelvic floor and slowly rotate your pelvis to the right, then to the left.

Repeat 20 times before lowering your pelvis to the floor and relaxing.

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
15-mal wiederholen.

## EN Pelvic push

Lie on your back, with your arms at your sides. Put the soles of your feet together and pull them up towards your crotch. Let your knees fall to the sides, and activate your pelvic floor. Contract your stomach muscles, pelvic floor, back, and buttocks, and slowly raise your pelvis. Relax your pelvic floor and lower your pelvis.

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## EN Waist kick

Lie on your left side, with your head resting on your extended left arm. Put a folded-up towel, blanket, or small cushion between your legs.. Place your right hand on the floor by your chest. Keep your torso and legs aligned and gently press your legs together. Activate your abdominal muscles, pelvic floor and back, and raise your legs from the floor. Raise your legs as high as you can. Lower them slowly without placing them back on the floor.

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
## EN Pelvic floor relaxation exercise

Roll up a blanket or keep a small cushion on hand, ready to use. Lie on your back and place the blanket or cushion in the small of your back. Let your feet fall to the sides in a relaxed position and stretch your arms out above your head. Hold this position for a few breaths and stretch the stomach and pelvic floor muscles before slowly relaxing them again.

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


## EN Bird dog position

Get down on your hands and knees. The knees should be directly beneath the hips, and the hands directly beneath the shoulders.

First, once again contract the pelvic floor. Then begin to slowly raise your right arm in front of you and your left leg behind you. Your arm, torso, and leg should form a straight line.

Make yourself as long as possible. Then slowly lower your arm and leg again before repeating the exercise on the other side.

Repeat a total of 10 to 12 times. Then rest before performing the exercise another 10 to 12 times.

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
## EN Bridge

Lie on your back, place your feet on the floor about hip-width apart, with your arms relaxed at your sides.

Tense your pelvic floor and slowly raise your pelvis and lower back off the floor. Your shoulder blades remain in contact with the floor. Try to keep your buttock muscles relaxed in the process, as this makes the pelvic floor work even harder. Hold this position for a few seconds and then gradually return to the starting position.

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Variation: Place your feet close together and press your knees together when raising the pelvis and back.

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## EN Prone position

Lie on your stomach and bend your right leg so that you are lying in a comfortable position. Place your forehead on top of your hands. How to carry out the exercise: Now pull in your stomach so that the pelvic floor contracts. Hold this position for a few seconds before once more relaxing.

Repeat a total of 8 times before changing legs, i.e. bend the left leg, and extend the right leg.

An intensive variation: Contract the pelvic floor, maintain the contraction, and then tense the pelvic floor yet further with 4 to 5 intense, quick contractions.

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
## EN Side lift

Turn onto your left side. Rest your left forearm on the floor so that your elbow is directly beneath your shoulder. Now initially push your shoulder away from your ear to “make yourself big”.

Now pull your navel in towards your spine, i.e. contract the pelvic floor. Maintain this tension and now slowly raise your hip off the floor.

Repeat the exercise about 6 times and then change sides.

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
## **EN *Knee raises***

Now raise your knee and press your right hand against the inside of your left knee as you exhale – the pelvic floor and the transverse abdominal muscles will contract. When you inhale, release the tension and repeat the exercise with your right leg.

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
## EN Cork screw

Lie on your back and stretch out both legs in a straight line towards the ceiling. Open your legs to hip-width, exhale, and place the left leg against the right leg. Keep your bottom in contact with the mat!

During this movement, your right foot is considerably lower than your left. Hold this position for a few long, deep breaths.

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## EN Side position

Lie on your side and draw your knees towards your stomach. Lightly contract your pelvic floor throughout the entire exercise.

Maintain your breathing flow and lift the upper leg so that your knee points towards the ceiling. Keep both feet together in the process.

Then bring the knees together from the hip and turn your foot towards the ceiling whilst keeping the knees together.

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
### EN Leg circles and heel beats

Lie on your stomach and place your forehead on your hands. Push your shoulders back, gently pull your tailbone towards your heels, causing the lower back to stretch. Hold this position throughout the entire exercise.

Slightly raise a stretched leg and rotate your heel inwards. Describe small circles with the leg (ten clockwise circles, ten anti-clockwise circles). In the same position, lift the other leg and repeat the exercise.

Back in the starting position, the back is stretched and the pelvic floor consciously contracted. Lift both legs slightly and gently “bang” your heels together ten times. Breathe in a continuous flow throughout.

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## EN Cross-legged sitting exercise

Sit down in a cross-legged position, keep your back straight, and prop yourself up with your hands. Now lift your pelvic floor muscles upwards and inwards.

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
## EN Plank exercise

Place your forearms and the tips of your toes on the floor. As you exhale, lift your legs and torso up off the mat – your pelvic floor will contract.

Hold this position for 15–20 seconds and repeat 3 times.

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
## EN *Kneeling exercise*

Kneel on the floor with your torso upright. Tense your body to maintain a good erect position (hand behind the head, elbows back, contracted lower back), with your upper (thoracic) spine upright, your sternum raised, and your neck stretched upwards. As you exhale, tilt your entire body slightly backwards (like a “tilting doll”) – your pelvic floor will contract.

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
## EN Relaxation exercise in prone position

Lie on the floor in a comfortable position, with a cushion under your pelvis and another under your head. Allow yourself to settle, and then relax completely. Now tense your arms by bending your elbows.

Exhale deeply and listen to the sound of your breathing. You can also make a gentle moaning or sighing noise. Now you can recognise and enjoy the difference between relaxation and tension.

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
## EN Chair exercise 1

Standing, with the legs about hip-width apart.

Contract your pelvic floor, maintain the tension, then cough before relaxing. Repeat this several times, and then try to lift the chair at the same time.

It is important to use correct lifting technique: Place your feet close to the chair, bend your knees and hips, then contract your pelvic floor, hold the tension, and lift the chair up by straightening your legs, before returning the chair to the floor. When you feel confident in using this lifting technique, lift the chair whilst simultaneously trying to cough.

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
## EN Chair exercise 2

Sit on a chair in an easy, relaxed position. With rounded shoulders, push your pelvis forwards. Exhale slowly whilst tightly squeezing shut your anal and vaginal passages. Maintain the tension for as long as you can. Count the seconds! Inhale, simultaneously releasing the tension. Your performance will improve each time you do this!

And pelvic floor training balls can bring success in less time and effort. Repeat the exercise 5 times for 5 seconds.

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
## EN Chair exercise 3

Sit down on a chair and inhale. Shift your weight in such a way that you are sitting on your vulva.

Exhale slowly and strongly contract the muscles of this orifice.

Count the seconds again: 10 seconds represents a very good performance. 15 would be excellent!

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
## EN Chair exercise 4

Sit upright on a chair and, as you inhale, stretch your arms out vertically above you. As you exhale, activate your pelvic floor. But also make sure you release the tension as you inhale. Slowly lower your arms. Now you can try out different arm positions. For example: hold them out in front of you, or to the left, or to the right.

Repeat the exercise 10 times for 5 seconds each time.

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

## EN Chair exercise 5

Sit on a chair, with your feet touching the floor. Now, with your knees together and your back straight, lift your feet off the floor. Exhale while doing this, and contract your buttock and abdominal muscles for a few seconds. As you inhale, place your feet back on the floor.

Repeat this exercise at least ten times.

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
## EN Chair exercise 6

Straddle a chair with a bolster pillow or something similar tucked between your legs. Now strongly contract the pelvic floor muscles around the bolster pillow. View the bolster pillow as something that is trying to escape and which you have to hang on to.

Maintain the tension as long as you can, and gradually try to increase the duration of the contraction each time. The rest you take in between contractions should last at least as long as the contractions themselves.

Do not hold your breath, but breathe in a steady rhythm!

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
## EN Chair exercise 7

Sit in a chair and lean forward slightly, with your weight towards the back of your pelvis. As you exhale, contract your anus, vagina, and urethra, and maintain the tension whilst continuing to breathe evenly.

Release the tension after 10–15 seconds. Repeat the exercise 8 times.

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## EN Chair exercise 8

Sit with an upright posture, stretch out your legs in front of you, and cross your ankles.

As you exhale, press the outer edges of your feet together whilst simultaneously contracting your pelvic floor and your buttock muscles. Maintain the tension for 10–15 seconds.

Continue to breathe steadily!

Repeat the exercise 6 times and take at least 5 seconds break between the exercises.

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### EN *Exercising routinely is the key to success*

It is recommended that you do the exercise(s) between three and ten times each day. As the exercises occur internally, i.e. they are not visible to others, you can in principle do them discreetly at any time, in any place. You can exercise at fixed times, or you can simply do the exercise(s) whenever you feel like it throughout the day. They will become as routine a part of daily life as going to bed or cleaning your teeth.

This way, your pelvic floor will remain in top condition, and bladder problems can be prevented or alleviated. Important information: The exercises listed and described here are not a substitute for visiting your doctor or any potentially necessary medical treatment. If in doubt, talk to your doctor and discuss with him whether and in what form these exercises are appropriate for you. Simply make a selection from this great variety of exercises, and create your own personal and individual training programme for varied yet challenging daily pelvic floor training.

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*recommended by midwives  
and gynaecologists*

MADE IN GERMANY

www.JOYDIVISION.de